# Proteomic analysis of Taenia solium excretionsecretion proteins

<u>B. Victor</u>, K. Kanobana, S. Gabriël, K. Polman, N. Deckers, P. Dorny, A. M. Deelder<sup>\*</sup> and M. Palmblad<sup>\*</sup>

> ITM (Antwerp, Belgium) \* LUMC (Leiden, The Netherlands)

### T. solium proteomics

- The project.
- Objectives of the study.
- Global methodology.
- What are the results so far?
- Discussion.
- What is next?



### PhD project

#### An FWO - ICONZ funded project

High throughput analysis of *Taenia* excretionsecretion proteins (ESPs) for improved understanding of host-parasite interactions and optimization of immunodiagnostic tools in *Taenia* solium cysticercosis.

### Objectives

2 major objectives:

I. Analyze the excretion-secretion proteome of *T*. solium to improve current knowledge of hostparasite interaction.

So far, no studies of the complete ES proteome.

#### Objectives

2. Optimization of immunodiagnostic tools in *T*. solium cysticercosis.

PROBLEM: antigen detection ELISA detects viable cysticerci, but is not species specific.

SOLUTION: find proteins unique for T. solium (or T. hydatigena) and make the ELISA more specific.

Starting point is the same for both objectives:
 A complete analysis of *T. solium* excretion-secretion proteins.

Step 2: Sample separation

Step I: Sample collection Step 3: Liquid chromatography & tandem mass spectrometry

Step 4: Data analysis

# Methodology Step I

- 5 pigs from Peru and Zambia.
- In vitro culture of cysticerci.
- Collect culture medium (= ESPs) after 24h and 48h.
- Concentrate ESPs (MWCO = IkDa).

Step 2: Sample separation

Step I: Sample collection Step 3: liquid chromatography & tandem mass spectrometry

Step 4: Data analysis

# Methodology Step 2

- Separate ESPs using I-D PAGE.
- Cut lane in 48 slices.
- Tryptic digest: proteins are cut into peptides.



Step 2: Sample separation

Step I: Sample collection Step 3: Liquid chromatography & tandem mass spectrometry

Step 4: Data analysis

# Methodology Step 3

- LUMC, Leiden.
- Separate peptides with liquid chromatography.
- On-line with 'ion trap' mass spectrometer.
- Result: dataset with experimental spectra.



Step 2: Sample separation

Step I: Sample collection Step 3: liquid chromatography & tandem mass spectrometry

Step 4: Data analysis

# Methodology Step 4

- Compare exp. spectra to spectra from a protein database with X!Tandem (in silico digest).
- Protein identification using validated tools (PeptideProphet, iProphet and ProteinProphet) to estimate accuracy of peptide/protein identifications.
- Ideal result: list of all proteins in the ESP sample.

#### However...

You need a protein database that contains all proteins likely to be in your sample.

T. solium protein database is incomplete.
 NCBI: ~ 270 proteins.

Genome is not (yet) known.

#### Solutions?

How to get around this limitation:
I. Sequence *T. solium* genome.
2. Supplement *T. solium* protein database with proteins from other helminths.

All Taenia (821) all Schistosoma (29,953), all Echinococcus (1,146), all Trichinella (16,325), Sus scrofa (1,388) and cRAP (112).

### What are the results so far?

Origin	# proteins
T. solium	27
Sequences from other helminths	32
S. scrofa	17

T. solium proteins: mostly proteins used in EITB.

Gene ontology classification: stress response, metabolism, detoxification, proteolysis, ...

#### Discussion

Unknown proteins in ESPs.
Unmatched spectra = unknown proteins?
Are helminth ES proteomes roughly the same?

Look to genome for differences?

B. Victor, K. Kanobana, S. Gabriël, K. Polman, N. Deckers, P. Dorny, A. M. Deelder and M. Palmblad. Proteomic analysis of *Taenia solium* metacestode excretion-secretion proteins. Accepted for publication in PROTEOMICS on March 5, 2012.

#### What is next?

 Search T. solium Expressed Sequence Tag (EST) library.

Transcriptomics: mRNA – cDNA – EST.

Find unassigned spectra and try de novo analysis.

Find peptide from spectrum.

Chance of finding 'new' peptides/proteins.

#### What is next?

T. solium genome project(s)?
 Derive proteins from DNA.

Analyze T. hydatigena ESPs and compare with T. solium (discover differential proteins).

2-D PAGE and look for differences.

### Thank you!

Vet. Helm. and Med. Helm. units (ITM).

- In Dr. Palmblad (LUMC, Leiden, The Netherlands).
- Oniversity of Zambia and the Universidad National Mayor de San Marcos (Lima, Peru).
- The Research Foundation Flanders (FWO) and the European Union's Seventh Framework Program (ICONZ) for financial support.